

DESI in-situ Calibration System status

Sonia Karkar (IN2P3/LPNHE) Laurent Le Guillou (Sorbonne Univ./LPNHE)

> Spectrograph Telecon Jun 5th, 2018



Talk outline

- Assemble
- Ship
- Unpack
- Fisrt check
- Interface with ICS
- Interface with the ring
- Ongoing at LPNHE
- Cooming up Next



Assembly in LPNHE









Prepare shipping at LPNHE





Prepare shipping at LPNHE





Trust DHL

- Send to NOAO
- NOAO will take it to Kitt Peak




Received at the Mayall (via NOAO)


















First check



First check

Video showing the 4 lamps turning ON one afer the other




Interface with ICS




- Klaus and Ann where at the Mayall to install the ICS and Co on the newly received computers
- This gave us the opportunity to try to control the calibartion boxes with the ICS

Interface with ICS: setup






Interface with ICS





Interface with ICS: HGAr ON






Interface with ICS: Kr ON





Interface with ICS: Ne and Xe





Interface with ICS

- Give their DESI ip adress to each PDU
 - Learn what to do if PDU network goes wrong: use consol serial link
- After some adjustement of ICS Control ON/OFF works
- Read back lamps state to display in the GUI works
- Read sensors works too

Interface with ICS


Homework at LPNHE:


3D print new black plastic covers with a hole for the serial console

















































Homework at LPNHE:
Produce black anodized side panels for the drawers
Because they show trough the (soon not) white grid



Ongoing and Coming up next:

- Cd spectral drawer : prototype done, 4 drawers to be assembled
- Continuum drawers :
 - Baseline : one drawer with halogen lamps and a blue filter to reduce the red part of the spectrum
 - Spectra taken on our spectro-photometric testbench, data analysis ongoing (uncertainties with the throughput of our monochromator).
 - Not much flux below 400 nm → add UV & blue LEDs to complete (tests on-going)

Continuum lamps: halogens + blue filter



- Baseline : one drawer with halogen lamps and a blue filter to reduce the red part of the spectrum
- Spectra taken on our spectrophotometric testbench, data analysis ongoing (uncertainties with the throughput of our monochromator).
- Not much flux below 400 nm; add UV & blue LEDs to complete (tests on-going)



Synthetic continuum with LEDs

